

Fitting & Wearing Your SenSlip


Please take the time and patience to achieve the best results from your new **SenSlip**. The package containing the **SenSlip** is reusable and can be used for storage.

To Start: You will need to be clean, dry and have some talcum powder to hand. If this is the first time you have used the **SenSlip** Undergarment, you may need to take a little extra time to make sure that you are getting the best and most comfortable results. When worn correctly the **SenSlip** should be of the utmost comfort. It is a tube constructed from pre-vulcanized latex, flocked with viscose rayon fibres. When fitted the tube becomes a double layer covering the glans. These layers are breathable.

On opening the package you will notice that the **SenSlip** is inside out showing the white flocked surface. This is correct. You will see that the **SenSlip** has an 'upside' or 'front' and an 'underside' or 'back'. To determine this there is a small bridge at the narrowest point molded into the product. This is on the underside or back of the **SenSlip** and will, when fitted, cover the underside of the glans. The upside or front has a slightly longer length at the ribbed end to cover the base of the flaccid penis.

With the use of talcum powder and with the white flocked fibres on the outside of the **SenSlip** gently stretch open the portion of the **SenSlip** which is to be fitted to the glans and pull it towards you to cover the entire glans. You will see that the main body of the **SenSlip** is now hanging forward of the penis. With the main body of the **SenSlip** hanging forward from the penis gently roll or pull back the body of the **SenSlip** so that the white flocked surface is in contact with the shaft of the penis. This may take several attempts but it will become second nature. The glans should now be covered with the double layer of the **SenSlip** Undergarment.

"take your time to get it right"



You may need to adjust the amount of 'overhang' or adjust the length of the **SenSlip** to fit the penis. It is recommended that you do not do this to start with but leave it for a while until you are more familiar with the feel of the product. When you are ready to make any adjustments to the length of the **SenSlip** you can follow the pattern at the end (ribbed end) with a sharp pair of scissors to achieve the correct length.

It is strongly recommended that you wear the **SenSlip** for just a few hours a day to start with. When you become more familiar with the product you will want to wear it all day. The **SenSlip** should be a 'snug' fit but never too tight. The **SenSlip** is designed to stay on the penis all day, you may find the need to make minor adjustments to it when visiting the bathroom.

Looking After Your SenSlip

To obtain the maximum benefit from your **SenSlip** Undergarment, please observe these important points:

- 1) The **SenSlip** Undergarment should only be worn during your conscious hours. It should not be worn whilst asleep.
- 2) The **SenSlip** is a very flexible undergarment. If an erection is experienced there will be enough expansion in the **SenSlip**. If for any reason discomfort is experienced at any time, the **SenSlip** should be removed.
- 3) When wearing the **SenSlip**, ensure that it does not get wet. If you are about to engage in any bathing or water sports, the **SenSlip** should be removed. If it does become wet or damp, remove it and allow it to dry. When dry, use talcum powder to keep it soft.
- 4) The **SenSlip** has been designed as an undergarment. It is not recommended that you use it as a sexual aid or toy. Damage to the product may occur.
- 5) Wash the **SenSlip** Undergarment by hand in warm water with a small amount of detergent. Rinse and then dry the **SenSlip** in a warm place overnight. When dry, use talcum powder to soften.
- 6) Be careful to not overstretch the **SenSlip** when fitting or removing. The life of the **SenSlip** can be greatly increased by treating it with care. Do not let the **SenSlip** come into contact with any oil based lubricants or ointments.
- 7) Do not use the **SenSlip** if you have any open flesh wounds or sores and never insert the **SenSlip** into any part of the body.