

Your-Skin Cone

User guide

Additional details at our "Instructions" page

TLC Tugger
Proudly helping men restore
www.TLCTugger.com

Ron Low

President 847 414-1692

TLCTugger@Juno.com



Step 1: Roll the skin of your flaccid penis forward to cover your glans



Step 2: When the glans is covered, bunch up your skin



Step 3: Wedge the bunched-up skin into the Your-skin Cone



Step 4: Check for uniformity. Here the exposed part of the glans is nicely centered



Step 5: Make sure the cone hangs on by itself. Wear it under briefs (not boxers) and to bed at night



Step 6: If the cone was invisible, this is how your skin would look with the cone properly fitted



Step 7: Be prepared to gain a healthy shine when your skin is pulled back

Step 8: Hygiene

- a) Wash the cone daily
- b) Remove the cone to urinate
- c) Wash your glans and inner skin only with clear water or a non-soap cleanser like Cetaphil liquid
- d) DO take before pictures
- e) After a month or two, take a jog without the cone. If your penis is irritated from banging around in your shorts, you'll know the cone is working.